

# FlexArm Dual Holder System



**FlexArm Dual Holder System  
provides unlimited positioning  
with static stabilization of  
HurriChem  
and laparoscopic camera**

# FlexArm Dual Holder System

## FlexArm Dual Holder System

FlexArm Single Arms (2),  
Swivel Ball Mount, and  
OR Table Rail Clamp  
**Order# PDT-5720**



## FlexArm Laparoscopic Camera Quick-Grip Tip

Safely and quickly secures  
endoscopes  
**Order# PDT-5705 - 5mm**  
**Order# PDT-5710 - 10mm**



## FlexArm HurriChem Instrument Clamp

Stainless Steel Instrument  
clamp attachment  
with Hex fitting  
**Order# PDT-5708**



Quick-Grip Tips and  
Instrument Clamp  
sold separately

### Specifications:

Post Height - 2 to 16" (5.1 - 40.6cm)  
Horizontal Reach - 19" (48.2cm)  
Flexible Reach - 16.5" (41.9cm)  
Tips: 3" (7.6cm) length (approx.)

Weight of complete set  
(without tray):  
10.4 lbs (4.72 kg)

## FlexArm Sterilization Tray **Order# PDT-5712**

